


# ORTHOPAEDICS




## IMPLANTS AND INSTRUMENTS



Siebvorschläge  
Set Suggestions  
Proposiciones de juegos  
Propositions pour boîtes  
Proposte di sets



**Grundsieb Knochen**  
**Basic Set Bone**  
**Juego básico para huesos**  
**Boîte de base à os**  
**Set di base per ossa**





**Grundsieb Knochen**  
**Basic Set Bone**  
**Juego básico para huesos**  
**Boîte de base à os**  
**Set di base per ossa**

**Part 1**

Item number	Description	Quantity
<b>01.22.10</b>	Scalpel blade · fig 10 · sterile · packing unit 100 pieces	1
<b>01.22.15</b>	Scalpel blade · fig 15 · sterile · packing unit 100 pieces	1
<b>01.27.03</b>	Scalpel handle · no 3 · standard · 12.5 cm	1
<b>01.27.13</b>	Scalpel handle · no 3L · long · 21 cm	1
<b>06.00.14</b>	Dressing forceps · 14.5 cm · standard	2
<b>06.00.20</b>	Dressing forceps · 20.0 cm · standard	2
<b>06.05.14</b>	Tissue forceps · 14.5 cm · 1 × 2 teeth · standard	2
<b>06.10.14</b>	Tissue forceps · 14.5 cm · 2 × 3 teeth	2
<b>06.21.12</b>	ADSON Tissue forceps · 12 cm · 1 × 2 teeth · smooth	2
<b>09.61.14</b>	REDON Guide needle · curved · 14 French	2
<b>09.83.08</b>	FRAZIER Suction cannula · 8 French · angled · 12.5 cm	1
<b>10.18.18</b>	MAYO-HEGAR Needle holder · 18 cm	2
<b>14.23.62</b>	UNIVERSAL HM Wire cutting scissors · 12 cm	1
<b>14.48.04</b>	Wire twisting forceps · 18 cm · round	1
<b>15.44.12</b>	HALSTED-MOSQUITO Forceps · 12.5 cm · straight	4
<b>15.45.12</b>	HALSTED-MOSQUITO Forceps · 12.5 cm · curved	4
<b>15.70.14</b>	KELLY Hemostatic forceps · 14 cm · straight	4
<b>15.71.14</b>	KELLY Hemostatic forceps · 14 cm · curved	4
<b>16.11.16</b>	ROCHESTER-PEAN Hemostatic forceps · 16 cm · curved	2
<b>16.12.16</b>	ROCHESTER-OCHSNER Hemostatic forceps · 1 × 2 teeth · 16 cm · straight	2
<b>16.12.20</b>	ROCHESTER-OCHSNER Hemostatic forceps · 1 × 2 teeth · 20 cm · straight	1
<b>17.55.13</b>	BACKHAUS Towel clamp · 13.0 cm	6
<b>17.57.14</b>	EDNA (LORNA) Towel clamp · 14 cm	2
<b>19.50.24</b>	FOERSTER Sponge forceps · 24 cm · straight	6
<b>22.34.01</b>	VOLKMANN Retractor · 1 tooth · sharp · 22 cm	1
<b>22.34.04</b>	VOLKMANN Retractor · 4 teeth · sharp · 18 mm · 22 cm	2
<b>22.35.04</b>	VOLKMANN Retractor · 4 teeth · blunt · 18 mm · 22 cm	2
<b>22.72.04</b>	ISRAEL Retractor · 4 teeth · blunt · 40 × 40 mm · 25 cm	2





**Grundsieb Knochen**  
**Basic Set Bone**  
**Juego básico para huesos**  
**Boîte de base à os**  
**Set di base per ossa**



23.63.80



26.29.07



26.60.35



26.64.26



26.64.84



26.08.77



27.95.16  
27.95.17



27.95.20



28.44.10  
28.44.15



28.45.20



28.61.08  
28.61.10



28.62.06  
28.62.10



28.64.20  
28.64.30



28.65.20



29.56.51  
29.56.54



30.38.69



31.07.05



33.20.17



35.11.17



35.14.16  
35.14.22



35.15.18



35.16.21



35.20.20



**Grundsieb Knochen**  
**Basic Set Bone**  
**Juego básico para huesos**  
**Boîte de base à os**  
**Set di base per ossa**

**Part 2**

Item number	Description	Quantity
<b>23.63.80</b>	HIBBS Retractor · toothed · set/2	1
<b>26.29.07</b>	LANGENBECK Elevator · 7 mm · 20 cm	1
<b>26.60.35</b>	LANGENBECK Raspatory · 15 mm · 19 cm	1
<b>26.64.26</b>	Raspatory · 6 mm · 18.5 cm	1
<b>26.64.84</b>	Raspatory · 14 mm · 20 cm · straight	1
<b>27.08.77</b>	OMBREDANNE Metal mallet · 770 g · 23 cm	1
<b>27.95.16</b>	CASPAR Pestle · Ø 5 mm · 20 cm	1
<b>27.95.17</b>	CASPAR Pestle · Ø 8 mm · 20 cm	1
<b>27.95.20</b>	CASPAR Pestle · Ø 12 mm · 20 cm	1
<b>28.44.10</b>	STILLE Osteotome · 10 mm · 20 cm	1
<b>28.44.15</b>	STILLE Osteotome · 15 mm · 20 cm	1
<b>28.45.20</b>	STILLE Gouge · 20 mm · 20 cm	1
<b>28.61.08</b>	MINI-LEXER Chisel · 8 mm · 18 cm	1
<b>28.61.10</b>	MINI-LEXER Chisel · 10 mm · 18 cm	1
<b>28.62.06</b>	MINI-LEXER Gouge · 6 mm · 18 cm	1
<b>28.62.10</b>	MINILEXER Gouge · 10 mm · 18 cm	1
<b>28.64.20</b>	LEXER Chisel · 20 mm · 22 cm	1
<b>28.64.30</b>	LEXER Chisel · 30 mm · 22 cm	1
<b>28.65.20</b>	LEXER Gouge · 20 mm · 22 cm	1
<b>29.56.51</b>	BRUNS Sharp curette · fig 1 · 27 cm	1
<b>29.56.54</b>	BRUNS Sharp curette · fig 4 · 27 cm	1
<b>30.38.69</b>	OLIVECRONA Bone rongeur · curved · 20 cm	1
<b>31.07.05</b>	RUSKIN Bone rongeur · 5 mm · curved · 18 cm	1
<b>33.20.17</b>	LISTON Bone cutting forceps · 17 cm · straight	1
<b>35.11.17</b>	KISTLER-MANUS Repositions forceps · 15.5 cm	1
<b>35.14.16</b>	Bone holding clamp · 16 cm	1
<b>35.14.22</b>	Bone holding clamp · 22 cm	1
<b>35.15.18</b>	DINGMAN Bone holding forceps · 19 cm	1
<b>35.16.21</b>	Malleolar holding clamp · 21 cm	1
<b>35.20.20</b>	Bone holding forceps · delicate · 20 cm · straight	1





**Grundsieb Knochen**  
**Basic Set Bone**  
**Juego básico para huesos**  
**Boîte de base à os**  
**Set di base per ossa**





**Grundsieb Knochen**  
**Basic Set Bone**  
**Juego básico para huesos**  
**Boîte de base à os**  
**Set di base per ossa**

**Part 3**

Item number	Description	Quantity
<b>35.45.78</b>	ULRICH Bone holding clamp · 28 cm · curved	1
<b>35.52.21</b>	FERGUSSON Bone holding clamp · 21 cm	1
<b>35.75.17</b>	VERBRUGGE Bone holding clamp · 17.5 cm	1
<b>35.75.25</b>	VERBRUGGE Bone holding clamp · 25 cm	1
<b>44.37.20</b>	Ruler · stainless steel · 20 cm	1
<b>47.20.25</b>	Kidney basin · 250 × 125 × 35 mm	1
<b>50.02.16</b>	ALLIS Grasping forceps · 15 cm · 5 × 6 teeth	4
<b>52.72.05</b>	SIMON Scoop · rigid · 29 cm · fig 5	1
<b>76.45.10</b>	Wire extraction pliers · serrated · 17 cm	1
<b>76.45.75</b>	GOLIATH HM Wire and plate cutting plier · 25 cm	1
<b>76.45.85</b>	EASY CUT HM Wire and pin cutter · 20 cm	1
<b>76.59.15</b>	Bone reamer · 15 cm · T handle	1
<b>77.10.10</b>	DINGMAN Bone hook · 10 mm · 19 cm	1





**Grundsieb Knochen**  
**Basic Set Bone**  
**Juego básico para huesos**  
**Boîte de base à os**  
**Set di base per ossa**



77.10.95



77.11.10  
77.11.18



77.11.24



77.11.42  
77.11.65



203.06.14



203.50.67



203.51.67



204.33.68



49.60.22



49.95.76



49.97.21  
49.97.22







**Grundsieb Knochen**  
**Basic Set Bone**  
**Juego básico para huesos**  
**Boîte de base à os**  
**Set di base per ossa**

**Part 4**

Item number	Description	Quantity
<b>77.10.95</b>	Bone lever · 25 mm · 27 cm	2
<b>77.11.10</b>	HOHMANN Bone lever · 10 mm · 22 cm	2
<b>77.11.18</b>	HOHMANN Bone lever · 18 mm · 24 cm	2
<b>77.11.24</b>	HOHMANN Bone lever · 24 mm · 26 cm	2
<b>77.11.42</b>	VERBRUGGE-HOHMANN Bone lever · 42 mm · 24 cm	2
<b>77.11.65</b>	VERBRUGGE-HOHMANN Bone lever · 65 mm · 25 cm	2
<b>203.06.14</b>	SC Operating scissors · straight · 14.5 cm · sharp / blunt	1
<b>203.50.67</b>	MAYO SC PLUS Scissors · straight · 17.0 cm	1
<b>203.51.67</b>	MAYO SC PLUS Scissors · curved · 17.0 cm	2
<b>204.33.68</b>	METZENBAUM SC PLUS Scissors · curved · 18 cm	1
<b>49.60.22</b>	Medicontainer · 590 × 290 × 220 mm · silver	1
<b>49.92.63</b>	Identification label red · engraved · 50 x 17.5 mm	1
<b>49.95.76</b>	Wrapping drape · 140 x 100 cm	1
<b>49.97.21</b>	Wire basket · 510 × 255 × 50 mm	1
<b>49.97.22</b>	Wire basket · 510 × 255 × 70 mm	2



**Grundsieb Trauma**  
**Basic Set Trauma**  
**Juego básico para trauma**  
**Boîte de base a trauma**  
**Set di base per trauma**





**Grundsieb Trauma**  
**Basic Set Trauma**  
**Juego básico para trauma**  
**Boîte de base a trauma**  
**Set di base per trauma**

**Part 1**

Item number	Description	Quantity
<b>01.22.10</b>	Scalpel blade · fig 10 · sterile · p/100	1
<b>01.22.15</b>	Scalpel blade · fig 15 · sterile · p/100	1
<b>01.27.03</b>	Scalpel handle · no 3 · standard · 12.5 cm	1
<b>01.27.13</b>	Scalpel handle · no 3L · long · 21 cm	1
<b>06.00.14</b>	Dressing forceps · 14.5 cm · standard	2
<b>06.00.20</b>	Dressing forceps · 20.0 cm · standard	2
<b>06.05.14</b>	Tissue forceps · 14.5 cm · 1 × 2 teeth · standard	2
<b>06.10.14</b>	Tissue forceps · 14.5 cm · 2 × 3 teeth	2
<b>06.21.12</b>	ADSON Tissue forceps · 12 cm · 1 × 2 teeth · smooth	2
<b>09.61.14</b>	REDON Guide needle · curved · 14 French	2
<b>09.83.08</b>	FRAZIER Suction cannula · 8 French · angled · 12.5 cm	1
<b>10.18.18</b>	MAYO-HEGAR Needle holder · 18 cm	2
<b>14.23.62</b>	UNIVERSAL HM Wire cutting scissors · 12 cm	1
<b>14.48.04</b>	Wire twisting forceps · 18 cm · round	1
<b>15.44.12</b>	HALSTED-MOSQUITO Forceps · 12.5 cm · straight	4
<b>15.45.12</b>	HALSTED-MOSQUITO Forceps · 12.5 cm · curved	4
<b>15.70.14</b>	KELLY Hemostatic forceps · 14 cm · straight	4
<b>15.71.14</b>	KELLY Hemostatic forceps · 14 cm · curved	4
<b>16.11.16</b>	ROCHESTER-PEAN Hemostatic forceps · 16 cm · curved	2
<b>16.12.16</b>	ROCHESTER-OCHSNER Hemostatic forceps · 1 × 2 teeth · 16 cm · straight	2
<b>16.12.20</b>	ROCHESTER-OCHSNER Hemostatic forceps · 1 × 2 teeth · 20 cm · straight	1
<b>17.55.13</b>	BACKHAUS Towel clamp · 13.0 cm	6
<b>17.57.14</b>	EDNA (LORNA) Towel clamp · 14 cm	2
<b>19.50.24</b>	FOERSTER Sponge forceps · 24 cm · straight	6
<b>22.34.01</b>	VOLKMANN Retractor · 1 tooth · sharp · 22 cm	1
<b>22.34.04</b>	VOLKMANN Retractor · 4 teeth · sharp · 18 mm · 22 cm	2
<b>22.35.04</b>	VOLKMANN Retractor · 4 teeth · blunt · 18 mm · 22 cm	2
<b>22.72.04</b>	ISRAEL Retractor · 4 teeth · blunt · 40 × 40 mm · 25 cm	2





**Grundsieb Trauma**  
**Basic Set Trauma**  
**Juego básico para trauma**  
**Boîte de base a trauma**  
**Set di base per trauma**





**Grundsieb Trauma**  
**Basic Set Trauma**  
**Juego básico para trauma**  
**Boîte de base a trauma**  
**Set di base per trauma**

**Part 2**

Item number	Description	Quantity
<b>17.55.10</b>	BACKHAUS Towel clamp · 10.5 cm	6
<b>19.40.26</b>	ULRICH Cotton swab forceps · 26 cm · straight	2
<b>20.12.75</b>	STRANDELL-STILLE Retractor · 17 cm · 6 × 15 mm	2
<b>22.35.04</b>	VOLKMANN Retractor · 4 teeth · blunt · 18 mm · 22 cm	2
<b>22.35.06</b>	VOLKMANN Retractor · 6 teeth · blunt · 28 mm · 22 cm	2
<b>22.35.08</b>	VOLKMANN Retractor · 8 teeth · blunt · 40 mm · 22 cm	2
<b>22.73.14</b>	BILLROTH Retractor · 4 teeth · 50 × 45 mm · 25 cm	2
<b>23.20.02</b>	KOCHER Retractor · 12 × 45 mm · 21 cm	2
<b>23.20.03</b>	KOCHER Retractor · 12 × 65 mm · 21 cm	2
<b>23.21.40</b>	KOCHER Retractor · 40 × 75 mm · 23 cm	2
<b>23.40.01</b>	KOCHER Retractor · 12 × 19 mm · 20 cm	2
<b>23.40.20</b>	KOCHER Retractor · 20 × 38 mm · 23 cm	2
<b>24.22.13</b>	WEITLANER Retractor · 3 × 4 teeth · sharp · 13.0 cm	1
<b>29.51.24</b>	VOLKMANN Sharp curette · 17 cm · fig 4/0	1
<b>29.51.33</b>	VOLKMANN Sharp curette · 17 cm · fig 3/0	1
<b>29.51.42</b>	VOLKMANN Sharp curette · 17 cm · fig 2/0	1
<b>29.51.50</b>	VOLKMANN Sharp curette · 17 cm · fig 0	1
<b>29.51.51</b>	VOLKMANN Sharp curette · 17 cm · fig 1	1
<b>29.51.52</b>	VOLKMANN Sharp curette · 17 cm · fig 2	1
<b>29.51.53</b>	VOLKMANN Sharp curette · 17 cm · fig 3	1
<b>29.51.54</b>	VOLKMANN Sharp curette · 17 cm · fig 4	1
<b>29.51.55</b>	VOLKMANN Sharp curette · 17 cm · fig 5	1
<b>29.51.56</b>	VOLKMANN Sharp curette · 17 cm · fig 6	1
<b>29.52.51</b>	SCHEDÉ Sharp curette · 17 cm · fig 1	1
<b>29.52.52</b>	SCHEDÉ Sharp curette · 17 cm · fig 2	1
<b>29.52.53</b>	SCHEDÉ Sharp curette · 17 cm · fig 3	1
<b>40.25.18</b>	LISTER Bandage scissors · 18 cm	1
<b>40.26.20</b>	LISTER Bandage scissors · 20 cm · large ring	1
<b>47.20.25</b>	Kidney basin · 250 × 125 × 35 mm	2
<b>47.31.08</b>	Round bowl · 150 × 75 mm · 1.00 l	2
<b>77.11.10</b>	HOHMANN Bone lever · 10 mm · 22 cm	2
<b>77.11.18</b>	HOHMANN Bone lever · 18 mm · 24 cm	4
<b>77.11.42</b>	VERBRUGGE-HOHMANN Bone lever · 42 mm · 24 cm	2
<b>78.35.11</b>	Bone lever · 6 mm · 16 cm	2





**Grundsieb Trauma**  
**Basic Set Trauma**  
**Juego básico para trauma**  
**Boîte de base a trauma**  
**Set di base per trauma**



87.33.70



87.32.04



87.32.20



87.32.28



88.00.21



88.00.40



88.20.08



204.33.18



49.60.22



49.95.76



49.97.21  
49.97.22




**Grundsieb Trauma**  
**Basic Set Trauma**  
**Juego básico para trauma**  
**Boîte de base a trauma**  
**Set di base per trauma**

**Part 3**

Item number	Description	Quantity
<b>87.32.04</b>	Spatula electrode · straight · Ø 4.0 mm	1
<b>87.32.20</b>	Wire loop electrode · Ø 5 mm · Ø 4.0 mm	1
<b>87.32.28</b>	Rollerball electrode · straight · Ø 4 mm · Ø 4.0 mm	1
<b>87.33.70</b>	Electrode handle · Ø 4.0 mm · with cable · 4.5 m	1
<b>88.00.21</b>	Bipolar cable · Ø 12.5 mm · 4 m	1
<b>88.00.40</b>	Bracket · for high frequency cables	1
<b>88.20.08</b>	Bipolar forceps · 20.0 cm · angled · blunt	2
<b>204.33.18</b>	METZENBAUM SC Scissors · curved · 18 cm	1
<b>49.60.22</b>	Medicontainer · 590 × 290 × 220 mm · silver	1
<b>49.92.63</b>	Identification label red · engraved · 50 x 17,5 mm	1
<b>49.95.76</b>	Wrapping drape · 140 x 100 cm	1
<b>49.97.21</b>	Wire basket · 510 × 255 × 50 mm	1
<b>49.97.22</b>	Wire basket · 510 × 255 × 70 mm	2



**Metallentfernung**  
**Removal of Plates and Screws**  
**Remoción de metal**  
**Extraction de métal**  
**Rimozione delle placche e viti**







**Metallentfernung**  
**Removal of Plates and Screws**  
**Remoción de metal**  
**Extraction de métal**  
**Rimozione delle placche e viti**

Item number	Description	Quantity
<b>01.22.10</b>	Scalpel blade · fig 10 · sterile · p/100	1
<b>01.22.22</b>	Scalpel blade · fig 22 · sterile · p/100	1
<b>01.27.03</b>	Scalpel handle · no 3 · standard · 12.5 cm	1
<b>01.28.04</b>	Scalpel handle · no 4 · standard · 13.5 cm	2
<b>06.05.14</b>	Tissue forceps · 14.5 cm · 1 × 2 teeth · standard	2
<b>09.60.12</b>	REDON Guide needle · angled · 12 French	1
<b>09.60.14</b>	REDON Guide needle · angled · 14 French	1
<b>09.85.00</b>	YANKAUER Suction tube · 28 cm	1
<b>19.17.26</b>	MAIER Dressing forceps · 26 cm · curved	2
<b>26.42.05</b>	WILLIGER Raspatory · 5 mm · 16 cm	1
<b>26.60.35</b>	LANGENBECK Raspatory · 15 mm · 19 cm	1
<b>31.07.25</b>	RUSKIN Bone rongeur · 5 mm · curved · 24 cm	1
<b>31.11.13</b>	MARQUARDT Bone rongeur · curved · 20 cm	1
<b>75.22.60</b>	Screwdriver set · 1 handle, 8 blades	1
<b>78.29.14</b>	BUNNELL Dissecting probe · 16.5cm · angled	1
<b>204.33.68</b>	SC Plus Metzenbaum Scissors · curved · 18 cm	1
<b>351.64.01</b>	Cleaning hook	1
<b>49.50.13</b>	Medicontainer · 290 × 290 × 135 mm · silver	1
<b>49.92.63</b>	Identification label red · engraved · 50 x 17.5 mm	1
<b>49.95.76</b>	Wrapping drape · 140 x 100 cm	1
<b>49.97.02</b>	Wire basket · 240 × 255 × 50 mm	2



**Ender Marknagelung**  
**Ender Intramedullary Nailing**  
**Ender fijación intramedular**  
**Ender enclouage centromédullaire**  
**Ender inchioldamento endomidollare**



75.66.30 - 78.66.49



75.66.05



75.66.06



75.66.07



75.66.08



75.95.21



75.95.61



49.60.15



49.97.21



49.95.76




**Ender Marknagelung**  
**Ender Intramedullary Nailing**  
**Ender fijación intramedular**  
**Ender enclouage centromédullaire**  
**Ender inchioldamento endomidollare**

Item number	Description	Quantity
<b>75.66.30</b>	ENDER Intramedullary hook · Ø 4.5 mm · 30 cm	3
<b>75.66.31</b>	ENDER Intramedullary hook · Ø 4.5 mm · 31 cm	3
<b>75.66.32</b>	ENDER Intramedullary hook · Ø 4.5 mm · 32 cm	3
<b>75.66.33</b>	ENDER Intramedullary hook · Ø 4.5 mm · 33 cm	3
<b>75.66.34</b>	ENDER Intramedullary hook · Ø 4.5 mm · 34cm	3
<b>75.66.35</b>	ENDER Intramedullary hook · Ø 4.5 mm · 35 cm	3
<b>75.66.36</b>	ENDER Intramedullary hook · Ø 4.5 mm · 36cm	3
<b>75.66.37</b>	ENDER Intramedullary hook · Ø 4.5 mm · 37 cm	3
<b>75.66.38</b>	ENDER Intramedullary hook · Ø 4.5 mm · 38 cm	3
<b>75.66.39</b>	ENDER Intramedullary hook · Ø 4.5 mm · 39 cm	3
<b>75.66.40</b>	ENDER Intramedullary hook · Ø 4.5 mm · 40 cm	3
<b>75.66.41</b>	ENDER Intramedullary hook · Ø 4.5 mm · 41 cm	3
<b>75.66.42</b>	ENDER Intramedullary hook · Ø 4.5 mm · 42 cm	3
<b>75.66.43</b>	ENDER Intramedullary hook · Ø 4.5 mm · 43 cm	3
<b>75.66.44</b>	ENDER Intramedullary hook · Ø 4.5 mm · 44 cm	3
<b>75.66.45</b>	ENDER Intramedullary hook · Ø 4.5 mm · 45 cm	3
<b>75.66.46</b>	ENDER Intramedullary hook · Ø 4.5 mm · 46 cm	3
<b>75.66.47</b>	ENDER Intramedullary hook · Ø 4.5 mm · 47 cm	3
<b>75.66.48</b>	ENDER Intramedullary hook · Ø 4.5 mm · 48 cm	3
<b>75.66.49</b>	ENDER Intramedullary hook · Ø 4.5 mm · 49 cm	3
<b>75.66.05</b>	Inserting Instrument for 75.66.30 - 49	1
<b>75.66.06</b>	Impactor for 75.66.30 - 49	1
<b>75.66.07</b>	Impactor/extractor for 75.66.30 - 49	1
<b>75.66.08</b>	Extracting Instrument for 75.66.30 - 49	1
<b>75.95.21</b>	KÜNTSCHER Awl · Ø 8 mm · 28 cm	1
<b>75.95.61</b>	Protecting shield for skin and tissue	1
<b>49.60.15</b>	Medicontainer · 590 × 290 × 150 mm · silver	1
<b>49.92.63</b>	Identification label red · engraved · 50 x 17.5 mm	1
<b>49.95.76</b>	Wrapping drape · 140 x 100 cm	1
<b>49.97.21</b>	Wire basket · 510 × 255 × 50 mm	2



**Basissieb Hüfte**  
**Basic Set Hips**  
**Set básico para cadera**  
**Boîte générale pour hanche**  
**Set di base per l'anca**





**Basissieb Hüfte**  
**Basic Set Hips**  
**Set básico para cadera**  
**Boîte générale pour hanche**  
**Set di base per l'anca**

**Part 1**

Item number	Description	Quantity
<b>01.22.22</b>	Scalpel blade · fig 22 · sterile · p/100	1
<b>01.22.24</b>	Scalpel blade · fig 24 · sterile · p/100	1
<b>01.28.04</b>	Scalpel handle · no 4 · standard · 13.5 cm	1
<b>01.28.14</b>	Scalpel handle · no 4L · long · 21 cm	1
<b>03.06.14</b>	Operating scissors · straight · 14.5 cm · sharp / blunt	1
<b>03.39.95</b>	RESANO Operating scissors · angled on flat · 25 cm	1
<b>04.33.70</b>	METZENBAUM HM Scissors · curved · 20 cm	1
<b>04.33.73</b>	NELSON-METZENBAUM HM Scissors · curved · 23 cm	1
<b>06.02.14</b>	Dressing forceps · 14.5 cm · slim	2
<b>06.05.16</b>	Tissue forceps · 16.0 cm · 1 × 2 teeth · standard	2
<b>06.05.20</b>	Tissue forceps · 20.0 cm · 1 × 2 teeth · standard	2
<b>06.06.20</b>	Tissue forceps · 20.0 cm · 1 × 2 teeth · medium	2
<b>06.64.19</b>	DEBAKEY AT Vascular forceps · 20 cm · 2.8 mm	2
<b>06.64.24</b>	DEBAKEY AT Vascular forceps · 24 cm · 2.8 mm	2
<b>09.60.12</b>	REDON Guide needle · angled · 12 French	1
<b>09.60.14</b>	REDON Guide needle · angled · 14 French	1
<b>09.60.16</b>	REDON Guide needle · angled · 16 French	1
<b>09.85.00</b>	YANKAUER Suction tube · 28 cm	1
<b>10.18.70</b>	MAYO-HEGAR HM Needle holder · 20 cm	2
<b>10.37.70</b>	BOZEMANN (WERTHEIM) HM Needle holder · 20 cm	1
<b>14.06.88</b>	HM Ligature scissors · 18.0 cm · curved	1
<b>15.41.12</b>	MICRO-HALSTED Hemostatic forceps · 12.5 cm · curved	6
<b>16.10.16</b>	ROCHESTER-PEAN Hemostatic forceps · 16 cm · straight	2
<b>16.12.16</b>	ROCHESTER-OCHSNER Hemostatic forceps · 1 × 2 teeth · 16 cm · straight	4
<b>16.12.22</b>	ROCHESTER-OCHSNER Hemostatic forceps · 1 × 2 teeth · 22 cm · straight	2
<b>16.41.18</b>	HALSTED (ADSON) Hemostatic forceps · 18 cm · curved	2





**Basissieb Hüfte**  
**Basic Set Hips**  
**Set básico para cadera**  
**Boîte générale pour hanche**  
**Set di base per l'anca**





**Basissieb Hüfte**  
**Basic Set Hips**  
**Set básico para cadera**  
**Boîte générale pour hanche**  
**Set di base per l'anca**

**Part 2**

Item number	Description	Quantity
<b>16.65.72</b>	OVERHOLT-GEISSENDOERFER Forceps · 22.5 cm	1
<b>17.55.10</b>	BACKHAUS Towel clamp · 10.5 cm	12
<b>19.16.26</b>	MAIER Dressing forceps · 26 cm · straight	2
<b>19.17.26</b>	MAIER Dressing forceps · 26 cm · curved	2
<b>22.16.16</b>	ROUX Retractor · 16.5 cm · fig 3	2
<b>22.43.06</b>	KOCHER Retractor · 6 teeth · blunt · 28 mm · 22 cm	2
<b>22.56.08</b>	KOERTE Retractor · 8 teeth · sharp · 24 cm	2
<b>22.72.06</b>	ISRAEL Retractor · 6 teeth · blunt · 60 × 50 mm · 25 cm	2
<b>23.25.12</b>	LANGENBECK-KOCHER Retractor · 20 × 60 mm · 21 cm	2
<b>23.69.05</b>	MICULICZ Retractor · 55 × 85 mm · 25.5 cm	2
<b>26.29.10</b>	LANGENBECK Elevator · 10 mm · 20 cm	1
<b>26.60.35</b>	LANGENBECK Raspatory · 15 mm · 19 cm	1
<b>26.67.34</b>	WAGNER Raspatory · 17 mm · 33 cm	1
<b>26.76.07</b>	SEMB Raspatory · 23 cm · fig 7	1
<b>27.10.70</b>	WAGNER Mallet · Ø 35 mm · 162/360 g · set complete	1
<b>27.95.17</b>	CASPAR Pestle · Ø 8 mm · 20 cm	1
<b>27.95.20</b>	CASPAR Pestle · Ø 12 mm · 20 cm	1





**Basissieb Hüfte**  
**Basic Set Hips**  
**Set básico para cadera**  
**Boîte générale pour hanche**  
**Set di base per l'anca**



28.44.60  
28.44.65

28.64.20  
28.64.25  
28.64.30

28.65.80  
28.65.85

28.65.93

28.94.30

29.57.53  
29.57.56

31.07.25

35.21.25

35.50.21



46.37.51

47.20.25

47.31.06  
47.31.08

76.59.22

77.10.20





**Basissieb Hüfte**  
**Basic Set Hips**  
**Set básico para cadera**  
**Boîte générale pour hanche**  
**Set di base per l'anca**

**Part 3**

Item number	Description	Quantity
<b>28.44.60</b>	STILLE Osteotome · curved · 10 mm · 20 cm	1
<b>28.44.65</b>	STILLE Osteotome · curved · 15 mm · 20 cm	1
<b>28.64.20</b>	LEXER Chisel · 20 mm · 22 cm	1
<b>28.64.25</b>	LEXER Chisel · 25 mm · 22 cm	1
<b>28.64.30</b>	LEXER Chisel · 30 mm · 22 cm	1
<b>28.65.80</b>	LEXER Gouge · 10 mm · curved · 25 cm	1
<b>28.65.85</b>	LEXER Gouge · 15 mm · curved · 25 cm	1
<b>28.65.93</b>	COBB Gouge · strong curved · 28 cm	1
<b>28.94.30</b>	LAMBOTTE Osteotome · straight · 30 mm · 24 cm	1
<b>29.57.53</b>	BRUNS Sharp curette · fig 3 · 33 cm	1
<b>29.57.56</b>	BRUNS Sharp curette · fig 6 · 33 cm	1
<b>31.07.25</b>	RUSKIN Bone rongeur · 5 mm · curved · 24 cm	1
<b>35.21.25</b>	BIRCHER-GANSKE Cartilage clamp · 20 cm	1
<b>35.50.21</b>	LANGENBECK Bone holding clamp · 20 cm	1
<b>46.37.51</b>	Inner case · only · for 46.37.50	1
<b>47.20.25</b>	Kidney basin · 250 × 125 × 35 mm	1
<b>47.31.06</b>	Round bowl · 110 × 55 mm · 0.40 l	1
<b>47.31.08</b>	Round bowl · 150 × 75 mm · 1.00 l	1
<b>76.59.22</b>	PERTHES Bone reamer · 22 cm · aluminum handle	1
<b>77.10.20</b>	DINGMAN Bone hook · 20 mm · 17 cm	1





**Basissieb Hüfte**  
**Basic Set Hips**  
**Set básico para cadera**  
**Boîte générale pour hanche**  
**Set di base per l'anca**



77.11.18



77.11.24  
77.11.34



77.11.42  
77.11.65



77.11.93



77.12.04



203.57.96



49.60.22



49.97.21  
49.97.22



49.95.76




**Basissieb Hüfte**  
**Basic Set Hips**  
**Set básico para cadera**  
**Boîte générale pour hanche**  
**Set di base per l'anca**

**Part 4**

Item number	Description	Quantity
<b>77.11.18</b>	HOHMANN Bone lever · 18 mm · 24 cm	2
<b>77.11.24</b>	HOHMANN Bone lever · 24 mm · 26 cm	2
<b>77.11.34</b>	HOHMANN Bone lever · 34 mm · 29 cm	2
<b>77.11.42</b>	VERBRUGGE-HOHMANN Bone lever · 42 mm · 24 cm	2
<b>77.11.65</b>	VERBRUGGE-HOHMANN Bone lever · 65 mm · 25 cm	2
<b>77.11.93</b>	MURPHY-LANE Bone lever · 31 cm	1
<b>77.12.04</b>	BENNETT Bone lever · 45 mm · 24 cm	1
<b>203.57.96</b>	LEXER SC PLUS Scissors · curved · 16 cm	1
<b>49.60.22</b>	Medicontainer · 590 × 290 × 220 mm · silver	1
<b>49.92.63</b>	Identification label red · engraved · 50 × 17,5 mm	1
<b>49.95.76</b>	Wrapping drape · 140 × 100 cm	1
<b>49.97.21</b>	Wire basket · 510 × 255 × 50 mm	1
<b>49.97.22</b>	Wire basket · 510 × 255 × 70 mm	2



**Hallux Valgus Korrektur**  
**Hallux Valgus Correction**  
**Corrección de los dedos de pie**  
**Correction d'orteil**  
**Correzione delle dita dei piedi**





**Hallux Valgus Korrektur**  
**Hallux Valgus Correction**  
**Corrección de los dedos de pie**  
**Correction d'orteil**  
**Correzione delle dita dei piedi**

**Part 1**

Item number	Description	Quantity
<b>01.22.10</b>	Scalpel blade · fig 10 · sterile · p/100	1
<b>01.22.22</b>	Scalpel blade · fig 22 · sterile · p/100	1
<b>01.27.03</b>	Scalpel handle · no 3 · standard · 12.5 cm	1
<b>01.28.04</b>	Scalpel handle · no 4 · standard · 13.5 cm	1
<b>03.57.96</b>	LEXER HM Scissors · curved · 16 cm	2
<b>04.33.68</b>	METZENBAUM HM Scissors · curved · 18 cm	1
<b>06.01.14</b>	Dressing forceps · 14.5 cm · medium	2
<b>06.05.14</b>	Tissue forceps · 14.5 cm · 1 × 2 teeth · standard	2
<b>06.07.14</b>	Tissue forceps · 14.5 cm · 1 × 2 teeth · slim	2
<b>09.60.10</b>	REDON Guide needle · angled · 10 French	2
<b>09.85.00</b>	YANKAUER Suction tube · 28 cm	1
<b>10.18.65</b>	MAYO-HEGAR HM Needle holder · 15 cm	1
<b>10.18.70</b>	MAYO-HEGAR HM Needle holder · 20 cm	1
<b>14.06.84</b>	HM Ligature scissors · 14.5 cm · curved	1
<b>14.48.05</b>	Seizing forceps · double action · 5 mm · 18 cm	1
<b>15.47.12</b>	HALSTED-MOSQUITO forceps · 12.5 cm · curved · 1 × 2 teeth	8
<b>16.10.14</b>	ROCHESTER-PEAN Hemostatic forceps · 14 cm · straight	4
<b>16.12.14</b>	ROCHESTER-OCHSNER Hemostatic forceps · 1 × 2 teeth · 14 cm · straight	6
<b>17.55.10</b>	BACKHAUS Towel clamp · 10.5 cm	2
<b>19.16.26</b>	MAIER Dressing forceps · 26 cm · straight	2
<b>19.17.26</b>	MAIER Dressing forceps · 26 cm · curved	2
<b>20.18.10</b>	LOVE (MORSE) Retractor · 10 mm · 25 cm	2
<b>22.16.14</b>	ROUX Retractor · 14.0 cm · fig 1	2
<b>22.30.02</b>	VOLKMANN Retractor · 2 teeth · sharp · 9 mm · 21 cm	2
<b>22.34.04</b>	VOLKMANN Retractor · 4 teeth · sharp · 18 mm · 22 cm	2





**Hallux Valgus Korrektur**  
**Hallux Valgus Correction**  
**Corrección de los dedos de pie**  
**Correction d'orteil**  
**Correzione delle dita dei piedi**



22.35.02 22.42.01 23.05.14 26.29.07 26.64.26P 27.08.57 27.82.14 28.61.08P 28.64.15P 28.65.10P 29.51.51



30.71.02 31.05.03 31.07.05 33.21.17 33.21.22 35.11.13 46.37.51 76.16.03



**Hallux Valgus Korrektur**  
**Hallux Valgus Correction**  
**Corrección de los dedos de pie**  
**Correction d'orteil**  
**Correzione delle dita dei piedi**

**Part 2**

Item number	Description	Quantity
<b>22.35.02</b>	VOLKMANN Retractor · 2 teeth · blunt · 9 mm · 22 cm	2
<b>22.42.01</b>	KOCHER Retractor · 1 tooth · sharp · 22 cm	1
<b>23.05.14</b>	LANGENBECK Retractor · 14 × 28 mm · 21 cm	2
<b>26.29.07</b>	LANGENBECK Elevator · 7 mm · 20 cm	1
<b>26.64.26</b>	Raspatory · 6 mm · 18.5 cm	1
<b>27.08.57</b>	BERGMANN (HIBBS) Metal mallet · 540 g	1
<b>27.82.14</b>	WEST Gouge · straight · 4 mm · 16 cm	1
<b>28.61.08</b>	MINI-LEXER Chisel · 8 mm · 18 cm	1
<b>28.64.15</b>	LEXER Chisel · 15 mm · 22 cm	1
<b>28.65.10</b>	LEXER Gouge · 10 mm · 22 cm	1
<b>29.51.51</b>	VOLKMANN Sharp curette · 17 cm · fig 1	4
<b>30.71.02</b>	BEYER Bone rongeur · 2 mm · curved · 17 cm	1
<b>31.05.03</b>	RUSKIN-MINI Bone rongeur · 3 mm · curved · 15 cm	1
<b>31.07.05</b>	RUSKIN Bone rongeur · 5 mm · curved · 18 cm	1
<b>33.21.17</b>	LISTON Bone cutting forceps · 17 cm · curved	1
<b>33.21.22</b>	LISTON Bone cutting forceps · 22 cm · curved	1
<b>35.11.13</b>	Bone clamp · 13.5 cm ang	2
<b>46.37.51</b>	Inner case · only · for 46.37.50	1
<b>76.16.03</b>	KIRSCHNER Wire tube · 16 cm · Ø 12 mm	1





**Hallux Valgus Korrektur**  
**Hallux Valgus Correction**  
**Corrección de los dedos de pie**  
**Correction d'orteil**  
**Correzione delle dita dei piedi**



76.45.85



77.10.20



77.10.60



77.11.10



78.35.11



87.30.12 1:1



87.33.70



204.33.63



49.60.17



49.95.76



49.97.21





**Hallux Valgus Korrektur**  
**Hallux Valgus Correction**  
**Corrección de los dedos de pie**  
**Correction d'orteil**  
**Correzione delle dita dei piedi**

**Part 3**

Item number	Description	Quantity
<b>76.45.85</b>	EASY CUT HM Wire and pin cutter · 20 cm	1
<b>77.10.20</b>	DINGMAN Bone hook · 20 mm · 17 cm	1
<b>77.10.60</b>	LANGE-BABY Bone lever · 17 mm · 17 cm	2
<b>77.11.10</b>	HOHMANN Bone lever · 10 mm · 22 cm	2
<b>78.35.11</b>	Bone lever · 6 mm · 16 cm	2
<b>87.30.12</b>	Lancet electrode · straight · 25 mm	1
<b>87.33.70</b>	Electrode handle · Ø 4.0 mm · with cable · 4.5 m	1
<b>204.33.63</b>	METZENBAUM SC PLUS Scissors · curved · 14 cm	1
<b>49.60.17</b>	Medicontainer · 590 × 290 × 170 mm · silver	1
<b>49.92.63</b>	Identification label red · engraved · 50 × 17,5 mm	1
<b>49.95.76</b>	Wrapping drape · 140 × 100 cm	1
<b>49.97.21</b>	Wire basket · 510 × 255 × 50 mm	3



**Orthopädie-Draht Set**  
**Orthopaedics Wire Set**  
**Set de alambres ortopédicos**  
**Jeu de fils orthopédique**  
**Set di fili ortopedico**



14.30.17



14.41.65



14.46.14



76.05.24  
76.05.27



76.16.04



76.45.19



76.45.21



76.45.30



76.45.73



76.01.10  
76.01.15



76.10.10  
76.10.15  
76.10.20  
76.10.25



49.60.10



49.97.12



**Orthopädie-Draht Set**  
**Orthopaedics Wire Set**  
**Set de alambres ortopédicos**  
**Jeu de fils orthopédique**  
**Set di fili orthopedico**

Item number	Description	Quantity
<b>14.30.17</b>	HM wire cutter · 17,5 cm	1
<b>14.41.65</b>	HM wire twisting forceps · 15 cm	2
<b>14.46.14</b>	Flat nose pliers · 14 cm	1
<b>76.05.24</b>	Wire guide · 24 cm	1
<b>76.05.27</b>	Wire guide · 27 cm	1
<b>76.16.04</b>	Wire container · Ø 12 mm · 19 cm	4
<b>76.45.19</b>	Flat nose pliers · 19 cm	1
<b>76.45.21</b>	Wire and Nail Extracting Forceps · 20 cm	1
<b>76.45.30</b>	Wire and Nail Extracting Forceps · 18 cm	1
<b>76.45.73</b>	HM wire cutter · 23 cm	1
<b>76.01.10</b>	Bone suturing wire · soft · Ø 1 mm · 10 m	1
<b>76.01.15</b>	Bone suturing wire · soft · Ø 1.5 mm · 10 m	1
<b>76.10.10</b>	Kirschner wire · Ø 1 mm · 15 cm · p/100	10
<b>76.10.15</b>	Kirschner wire · Ø 1.5 mm · 15 cm · p/100	20
<b>76.10.20</b>	Kirschner wire · Ø 2.0 mm · 15 cm · p/100	10
<b>76.10.25</b>	Kirschner wire · Ø 2.5 mm · 15 cm · p/100	10
<b>76.16.04</b>	Wire container · Ø 12 mm · 19 cm	4
<b>49.60.10</b>	Medicontainer · 590 x 290 x 100 mm · silver	1
<b>49.92.63</b>	Identification label 50 x 17.5 mm · engraved · red	1
<b>49.97.12</b>	Wire basket · 480 x 255 x 50 mm	1